


EXPERIMENT 1

AIM: To verify Lambert – beer 's law for KMnO_4 colorimetrically.

THEORY

The primary objective of this experiment is to determine the concentration of an unknown KMnO_4 solution. The KMnO_4 solution used in this experiment has a blue color, so Colorimeter users will be instructed to use the red LED. A higher concentration of the colored solution absorbs more light (and transmits less) than a solution of lower concentration.



You will prepare five of known concentration (standard solutions).

Each solution is transferred to a small, rectangular cuvette that is placed into the Colorimeter or Spectrometer. The amount of light that penetrates the solution and strikes the photocell is used to compute the absorbance of each solution. When you graph absorbance vs. concentration for the standard solutions, a direct relationship should result. The direct relationship between absorbance and concentration for a solution is known as *Beer's law*.

You will determine the concentration of an unknown KMnO_4 solution by measuring its absorbance. By locating the absorbance of the unknown on the vertical axis of the graph, the corresponding concentration can be found on the horizontal axis. The concentration of the unknown can also be found using the slope of the Beer's law curve.

REQUIREMENTS

Colorimeter cuvette

five 20×150 mm test tubes

two 10 mL pipets or graduated cylinders two 100 mL beakers

0.01M KMnO_4 solution

distilled water test tube rack stirring rod

tissues (preferably lint-free)

PROCEDURE

1. Obtain small volumes of 0.01M KMnO_4 solution and distilled water in separate beakers.
2. Label five clean, dry, test tubes 1–5. Use pipets to prepare five standard solutions according to the chart below. Thoroughly mix each solution with a stirring rod. Clean and dry the stirring rod between uses.

| Test Tube | 0.01M KMnO_4 (mL) | Distilled H_2O (mL) | Concentration (M) |
|-----------|----------------------------|-------------------------------------|-------------------|
| 1 | 2 | 8 | 0.002 |
| 2 | 4 | 6 | 0.004 |
| 3 | 6 | 4 | 0.006 |
| 4 | 8 | 2 | 0.008 |
| 5 | ~10 | 0 | 0.0100 |

4. Prepare a *blank* by filling a cuvette 3/4 full with distilled water. To correctly use cuvettes, remember:
 - Wipe the outside of each cuvette with a lint-free tissue.
 - Handle cuvettes only by the top edge of the ribbed sides.
 - Dislodge any bubbles by gently tapping the cuvette on a hard surface.
 - Always position the cuvette so the light passes through the clear sides.

You are now ready to collect absorbance-concentration data for the five standard solutions.

- a. Using the solution in Test Tube 1, rinse the cuvette twice with ~1 mL amounts and then fill it 3/4 full. Wipe the outside with a tissue and place it in the device
- b. ~~(Colorimeter or Spectrometer). Close the lid on the Colorimeter. Test Tube 2, (close and~~
- d. Repeat the procedure for Test Tubes 3 and 4. Trial 5 is the original 0.01M KMnO_4 solution. **Note:** Do not test the unknown solution until Step 9.
- e. When you have finished testing the standard solutions
- f.

-
10. Determine the absorbance value of the unknown KMnO_4 solution.
- Obtain about 5 mL of the *unknown* KMnO_4 in another clean, dry, test tube. Record the number of the unknown in your data table.
 - Place the test tube vertically into the device. (Close the lid of the Colorimeter.)
 - Read the absorbance value displayed in the meter. When the displayed absorbance value stabilizes, record its value as Trial 6 in your data table.
 - Select Interpolate from the Analyze menu. Find the absorbance value that is closest to the absorbance reading you obtained in Step c above. Determine the concentration of your unknown KMnO_4 solution and record the concentration in your data table.
 - Dispose of any of the remaining solutions as directed.

DATA TABLE


| Trial | Concentration (mol/L) | Absorbance |
|-------|-----------------------|------------|
| 1 | 0.002 | |
| 2 | 0.004 | |
| 3 | 0.006 | |
| 4 | 0.008 | |
| 5 | 0.010 | |
| 6 | Unknown number ____ | |

EXPERIMENT 4

AIM: To verify Lambert – Beer's law for $K_2Cr_2O_7$ colorimetrically.

THEORY

The primary objective of this experiment is to determine the concentration of an unknown $K_2Cr_2O_7$ solution. The $K_2Cr_2O_7$ solution used in this experiment has a blue color, so Colorimeter users will be instructed to use the red LED. A higher concentration of the colored solution absorbs more light (and transmits less) than a solution of lower concentration.



You will prepare five of known concentration (standard solutions).

Each solution is transferred to a small, rectangular cuvette that is placed into the Colorimeter or Spectrometer. The amount of light that penetrates the solution and strikes the photocell is used to compute the absorbance of each solution. When you graph absorbance vs. concentration for the standard solutions, a direct relationship should result. The direct relationship between absorbance

and concentration for a solution is known as *Beer's law*.

You will determine the concentration of an unknown $K_2Cr_2O_7$ solution by measuring its absorbance. By locating the absorbance of the unknown on the vertical axis of the graph, the

corresponding concentration can be found on the horizontal axis. The concentration of the unknown can also be found using the slope of the Beer's law curve.

REQUIREMENTS

Colorimeter cuvette

five 20 × 150 mm test tubes

two 10 mL pipets or graduated cylinders two 100 mL beakers

0.01M $K_2Cr_2O_7$ solution

distilled water test tube rack stirring rod

tissues (preferably lint-free)

DAV COLLEGE JALANDHAR

PROCEDURE

1. Obtain small volumes of 0.01M $K_2Cr_2O_7$ solution and distilled water in separate beakers.
2. Label five clean, dry, test tubes 1–5. Use pipets to prepare five standard solutions according to the chart below. Thoroughly mix each solution with a stirring rod. Clean and dry the stirring rod between uses.

| Test Tube | 0.01M $K_2Cr_2O_7$ (mL) | Distilled H ₂ O (mL) | Concentration (M) |
|-----------|-------------------------|---------------------------------|-------------------|
| 1 | 2 | 8 | 0.002 |
| 2 | 4 | 6 | 0.004 |
| 3 | 6 | 4 | 0.006 |
| 4 | 8 | 2 | 0.008 |
| 5 | ~10 | 0 | 0.0100 |

4. Prepare a *blank* by filling a cuvette 3/4 full with distilled water. To correctly use cuvettes, remember:
 - Wipe the outside of each cuvette with a lint-free tissue.
 - Handle cuvettes only by the top edge of the ribbed sides.
 - Dislodge any bubbles by gently tapping the cuvette on a hard surface.
 - Always position the cuvette so the light passes through the clear sides.

You are now ready to collect absorbance-concentration data for the five standard solutions.

- a. Using the solution in Test Tube 1, rinse the cuvette twice with ~1 mL amounts and then fill it 3/4 full. Wipe the outside with a tissue and place it in the device (Colorimeter or Spectrometer). Close the lid on the
- b. ~~Colorimeter. Use the cuvette in Test Tube 2. Wipe the outside of the cuvette.~~
- d. Repeat the procedure for Test Tubes 3 and 4. Trial 5 is the original 0.01M $K_2Cr_2O_7$ solution. **Note:** Do not test the unknown solution until Step 9.
- e. When you have finished testing the standard solutions
- f.

DAV COLLEGE JALANDHAR

10. Determine the absorbance value of the unknown $K_2Cr_2O_7$ solution.

- a. Obtain about 5 mL of the *unknown* $K_2Cr_2O_7$ in another clean, dry, test tube. Record the number of the unknown in your data table.
- b. ~~Rinse the cuvette with the unknown solution and fill it (Close the~~
Rinse the cuvette with the unknown solution and fill it (Close the
- c. Read the absorbance value displayed in the meter. When the displayed absorbance value stabilizes, record its value as Trial 6 in your data table.
- d. Select Interpolate from the Analyze menu. Find the absorbance value that is closest to the absorbance reading you obtained in Step c above. Determine the concentration of your unknown $K_2Cr_2O_7$ solution and record the concentration in your data table.
- e. Dispose of any of the remaining solutions as directed.

DATA TABLE

| Trial | Concentration (mol/L) | Absorbance |
|-------|-----------------------|------------|
| 1 | 0.002 | |
| 2 | 0.004 | |
| 3 | 0.006 | |
| 4 | 0.008 | |
| 5 | 0.010 | |

| | | |
|---|-------------------------|--|
| 6 | Unknown number _____ | |
|---|-------------------------|--|

DAV COLLEGE JALANDHAR

DAV COLLEGE JALANDHAR